

COLLOÏDS & INTERFACES

**CEMENT, MORTAR & CONCRETE** 

# **ZETACAD:** ZETAMETRY USING STREAMING POTENTIAL OR CURRENT



#### **Parameters Measured**

- Streaming Potential
- > Streaming current (optional)
- > Plug resistance
- > Electrical Conductivity
- > Temperature

#### **Features and Benefits**

- > Applicable to particles above 50 μm diameter and flat surfaces.
- Reliable and simple to setup. Measurement and rinsing of the system are fully automated.
- Menu driven software Windows based.
- Data acquisition creates ASCII files which are directly compatible with common spreadsheets.


# ZETA POTENTIAL MEASUREMENT USING THE POROUS PLUG TECHNIQUE

An electrolyte is forced to pass through a capillary or porous plug by a pressure gradient. The excess charges around the particles or wall are carried along by the liquid.

Their accumulation downstream causes the build-up of an electric field which drives an electric current back (by ionic conduction) through the liquid, against the direction of the liquid flow.

A steady state is quickly established and the measured potential across the capillary is called the streaming potential.

#### Streaming current/potential of Fluorine





# **ZETACAD: ZETAMETRY USING STREAMING POTENTIAL OR CURRENT**

CAD Instruments offers wide range of services to help you take advantage of this new measurement device. The **ZetaCAD** can be used for major industrial and academic applications including:

- > Ceramics
- > Fibre and Textile
- > Membranes

- > Water treatment
- > Pulp & Paper

- > Polymers
- > Geology

#### ZETACAD® SPECIFICATIONS

## **Technology**

> Zeta Potential analysis by streaming potential determination using the porous plug technique

### **Specifications Range**

> Differential Pressure	± 500	mbar	
> Streaming Potential	± 2400	mV	
> Streaming Current	± 240	μΑ	

#### **Measuring Cell**

> Standard cell diameter (others upon request)	15	mm
> Variable cell length	10 -150	mm
> Tangential cell (surface specimen chamber	40 mm x 50 mm x 10 mm (W x D x H)	
> Flow through cell diameter (variable thickness)	47	mm

#### **Characteristics**

> Conductivity meter	0 – 20	mS.cm <sup>-1</sup>	
> Power Supply	100 – 240 V // 50 – 60 Hz // 50 VA		
> Dimensions	600 mm x 600 mm x 600 mm (W x D x H)		
> Weight	40	kg	
> Communication	USB bi-directional interface		

## Requirements

> Minimum computer specification: Pentium 4 with Windows® 2000, XP or Seven – 512 MB

Note: These specifications may change in the interest of product development The ZetaCAD was designed in cooperation with University of Nancy, France

